Instructions for use

Sophy® Mini Monopressure Valve

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Instructions for use

SOPHY[®] MINI MONOPRESSURE VALVE

Fixed pressure valve for CSF shunting

Sterile, single use

EN

Table of content

1.	Indications	5
2.	Contraindications	5
3.	Description and Operating Principle of the Sophy® Mini monopressure valve	5
4.	Configurations of the Sophy® Mini monopressure valve	5
5.	Measurement Unit and Calibration of Operating Pressures	6
6.	Sterilisation – Decontamination of Valves and Sophy® Mini monopressure valve Kits	6
7.	Instructions	6
	7.1. Choice of valve model	6
	7.2. Implantation Technique	6
	7.3. Post-operative X-ray control: Identification of the valve model and pressure reading	7
	7.4. Patency Test (post-operative)	8
	7.5. Sampling the CSF and injection	8
8.	Precautions for the Daily Life of the Patient	8
9.	Complications / Side effects	8
	9.1. Obstruction	9
	9.2. Infection	9
	9.3. Overdrainage	9
	9.4. Other	9
10.	Storage conditions	10
11.	Processing of the products after use	10
	11.1. Destruction after use	10
	11.2. Return of products	10
12.	Warranty	10
13.	Symbols	10
14.	References	10
	14.1. Sophy® Mini monopressure valves for CSF shunting	10
	14.2. Shunt kits for cranial implantation	11

All the illustrations are given for a medium pressure valve.

The same configurations exist in low and high pressure.

FIGURE 1 – SOPHY® MINI MONOPRESSURE VALVE (SM1, SM1A, SM1B). TOP AND SIDE VIEWS.



FIGURE 2 – COMPOSITION OF A COMPLETE SOPHY® Mini SM1-2010 (WITHOUT RESERVOIR) AND SM1A-2010 (WITH ANTECHAMBER RESERVOIR) KIT.

WARNING

Federal (USA) Law restricts this device to sale by or on the order of a physician.

1. Indications

The Sophy® Mini monopressure valve is intended for the treatment of hydrocephalus by shunting the Cerebrospinal Fluid (CSF) to the abdominal cavity or the right atrium of the heart.

2. Contraindications

The Sophy® Mini monopressure valves are not designed, sold or offered for use for any purpose other than their indication.

The contra-indications are the following :

- established or suspected infections along the length of the shunt (meningitis, ventriculitis, peritonitis, septicemia or bacteremia) or any infection present in any part of the body whatsoever,
- patients on anticoagulant therapy, or presenting with bleeding diathesis, or with hemorrhagic CSF, as the presence of blood in the system could lead to an obstruction in the system),
- ventriculo-atrial shunts in patients suffering from congenital cardiopathies or other malformations of the cardio-pulmonary system.

WARNING

The use of an external shunting devise (collection bag, etc.) in series with a valve should be avoided. However, in the case of temporary externalisation of a ventriculo-atrial/peritoneal shunt fitted with a SOPHY® Mini Monopressure valve, monitoring must be carried out by competent personnel in accordance with the procedures normally used for external ventricular shunts.

Although possible, the abovementioned cases do not represent normal use. If the clinical benefit is considered highly significant, their implementation is carried out under the responsibility of an informed neurosurgeon. The clinical condition of the patient must be subject to increased surveillance.

3. Description and Operating Principle of the Sophy® Mini monopressure valve


FIGURE 3: THE SOPHY® MINI MONOPRESSURE VALVE

The Sophy® Mini monopressure valve (Figures 1 & 3) is a single use implantable device.

The Sophy® Mini monopressure valve makes unidirectional drainage of the Cerebrospinal Fluid (CSF) possible.

The CSF arrives in the valve through the inlet connector [1], passes through the body of the valve [3] and leaves through the outlet connector [8].

The connectors are made of stainless steel and the body of the valve is in polysulfone.

On either side of the inlet connector [1] are suture holes [5] which make it possible to attach the valve to the subcutaneous tissues to prevent the valve from migrating.

On the upper surface of the valve an arrow shows the direction of CSF flow through the valve. This helps to position the valve correctly during implantation.

On the lower surface of the valve there is an individual serial number.

The valve body contains a ball-in-cone mechanism which determines the operating pressure of the valve.

The principle is based on the pressure exerted on a ruby ball [2] by a flat, semi-circular spring [6].

In conditions of normal use this mechanism provides an anti-reflux function and is not sensitive to temperature variations.

The valve body, which cannot be deformed, protects the mechanism from mechanical shocks.

It also prevents attempts at pumping or puncturing the valve and makes it insensitive to variations in percutaneous pressure.

The Sophy® Mini monopressure valve does not contain phthalates, or natural or synthetic latex.

Titanium radio-opaque studs [7] inserted on the right-hand side of the valve body indicate the operating pressure value for a given fixed pressure Sophy® Mini valve (See § - Post-operative X-ray control).

4. Configurations of the Sophy® Mini monopressure valve

There are 3 different pressures for the Sophy® Mini monopressure valve, Low (L), Medium (M) and High (H).

Sophy® Mini monopressure valves are available in 3 models:

- without a reservoir (SM1),
- with an integral reservoir of the antechamber type (SM1A),
- with an integral reservoir of the burr hole type (SM1B).

Models SM1 and SM1A are available as the valve alone or as a complete kit (Figures 1 & 2).

These configurations are detailed in the following table.

Operating pressures

References concerned	SM1-L SM1-2010L SM1A-L SM1A-2010L SM1B-L	SM1-M SM1-2010M SM1A-M SM1A-2010M SM1B-M	SM1-H SM1-2010H SM1A-H SM1A-2010H SM1B-H
Pressure (mmH ₂ O)	50 Low	110 Medium	170 High

Sophysa offers a complete range of radio-opaque ventricular and distal catheters which allow the CSF to flow to the valve and from the valve to the peritoneum or right atrium respectively, depending upon the type of shunt chosen by the neurosurgeon.

To be complete a Sophy® Mini monopressure shunt system must consist of a ventricular catheter, a Sophy® Mini valve and a distal catheter (atrial or peritoneal).

5. Measurement Unit and Calibration of Operating Pressures

The pressures mentioned are in mmH₂O.

1 mmH₂O corresponds to 9.807 Pa or 0.074 mmHg.

The valves are calibrated on the basis of a flow rate of 10 ml/h.

Each valve is tested individually: the measurement concerns the upstream pressure of a 10 ml/h flow of water passing through the valve and the Sophysa proximal and distal catheters.

The calibration is performed disregarding the resistance of the catheters.


Thus the pressures given on the valve labels correspond to the resistance of the valve alone.

The catheters add their own resistance to the shunt.

The reservoirs are not considered to have any particular resistance.

The calibration of Sophy® Mini valves is performed with a tolerance of -10/+15 mmH₂O on the measured pressures.

FIGURE 4: FLOW RATE-PRESSURE CURVES FOR THE SM1 MODEL

This curve is obtained for each pressure level (H, M, L) by varying the applied pressure and measuring the flow rate obtained. The values are given disregarding the resistance of the catheters.

Table of mean values at 10ml/h for each position setting of the Sophy® Mini valve, and corresponding values at 20ml/h:

	at 10 ml/h (mmH ₂ O)	at 20 ml/h (mmH ₂ O)
SM1-L	50	55
SM1-M	110	125
SM1-H	170	190

The effects on the device operating pressure of changes in the position of the patient and of sub-cutaneous pressure are negligible.

6. Sterilisation – Decontamination of Valves and Sophy® Mini monopressure valve Kits

The valves and Sophy® Mini monopressure valve kits are packed individually in double peel-off, sterile, pyrogen-free packaging, and have been sterilised with ethylene oxide.

Do not use the valves or valve kits if the sterile packaging is open or damaged, or if the expiry date has passed.

This product is intended for SINGLE USE ONLY. It is intended to be used once only for a single patient. Do not re-sterilize or re-use after unpacking and/or explantation.

Resterilization can damage the product, potentially leading to patient injury. Reuse of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

NOTE

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized and/or re-used, or for any complications which might result from this.

7. Instructions

7.1. CHOICE OF VALVE MODEL

The choice of Sophy® Mini monopressure valve model is left to the initiative of the neurosurgeon, depending on the clinical needs of the patient.

NOTE

Use an SM1A (antechamber), an SM1B (burr hole reservoir) model, or one of the SM1 models combined with a ventricular catheter with reservoir if the use of the shunt system is desired to check the patency of the shunt, to sample the CSF or for injections.

7.2. IMPLANTATION TECHNIQUE

Implantation of a Sophy® Mini monopressure valve must take account of current aseptic neurosurgical practices.

The implantation of a shunt including a Sophy® Mini monopressure valve may be performed in several ways.

The surgeon will choose the technique depending upon his experience and the clinical status of the patient.

The final implantation of the device must satisfy the conditions for optimal drainage of the CSF.

NOTE

Do not carry out an additional test before implantation: each valve has been individually calibrated and checked. Any pre-operative pressure tests will increase the risk of infection.

PRECAUTION

Avoid any contact of devices with contaminating items.

7.2.1. Ventricular Catheter

- Introduce the catheter into the ventricle using the introducing stylet supplied for this purpose.
- If necessary, adjust the implantation depth of the ventricular catheter with the right angle adapter supplied. Position it in the axis of the burr hole.
- Purge the catheter of air with the CSF.
- If necessary, check that the reservoir is properly filled, and then clamp.
- Connect and delicately ligate the catheter to the inlet connector of the valve (or that of the reservoir for valve models with integrated reservoirs). Check that the arrow located on the upper surface of the valve is correctly oriented in the direction of the flow. The clamp can then be released.

PRECAUTION

Ensure that the arrow on the upper surface of the valve is correctly oriented in the direction of the flow: assembly of the valve in the opposite direction would prevent any drainage.

PRECAUTION

Due to the fragility of the silicone, using metallic forceps for inserting catheters and for ligating them onto connectors is not recommended. This would create a risk of cutting or piercing the catheters.

7.2.2. Valve

WARNING

In the case of valves with an integrated reservoir or catheter, do not attempt to detach the reservoir or catheter from the valve. Detachment of the reservoir or catheter may unscrew the connector closure screw and uncalibrate the valve.

Sophy® Mini monopressure valves may be implanted either on the skull or in the supraclavicular fossa. The latter implantation may be preferred in paediatrics to avoid any risk of skin erosion.

1. Purge the valve of air. To prevent any risk of introducing an air bubble, it is recommended that the valve be left to fill directly with the patient's CSF. In the majority of cases, the valve fills immediately.

However, in patients with low intracranial pressure or if the valve has a high pressure value, the valve may not fill spontaneously. In this case :

- place a piece of catheter on the outlet connector and slowly aspirate the CSF using a syringe fitted with a Luer connector,
or even:
- press the dome of the reservoir to allow the CSF to fill the valve (SM1A and SM1B models).

PRECAUTION

Do not fill or purge the valve with any liquid other than the patient's CSF or water for injection (WFI) before implantation to avoid any risk of deposits in the valve, which could lead to an obstruction in the shunt system or a blockage in the valve mechanism.

2. Check that the valve is correctly filled with CSF and there are no air bubbles inside the valve. If this is not the case, continue to purge. The presence of air bubbles could cause a significant change to the operating pressure initially chosen.
3. Check that the arrow on the upper surface of the valve is visible and correctly oriented in the direction of the CSF flow.

NOTE

Do not implant the valve without suturing it to the underlying tissues by its two connectors or by the suture holes provided for this. If the shunt system migrates the drainage may stop and other complications ensue.

7.2.3. Peritoneal Catheter

1. Make a short peritoneal incision in the peri-umbilical region.
2. Tunnel the distal catheter.
3. Connect the proximal end of the catheter to the valve outlet connector and ligate it delicately.
4. Check the flow of CSF.
5. Adapt the length of the catheter.
6. Bury the distal end of the catheter in the peritoneal cavity.

PRECAUTION

Due to the fragility of the silicone, using metallic forceps for inserting catheters and for ligating them onto connectors is not recommended. This would create a risk of cutting or piercing the catheters.

7.3. POST-OPERATIVE X-RAY CONTROL: IDENTIFICATION OF THE VALVE MODEL AND PRESSURE READING




NOTE

During the radiological examination orientate the patient so that the ray source points perpendicularly onto the valve body.

In this way identifying the valve by its radio-opaque point is made easy.

Locate the valve inlet connector, wider due to the presence of a nut.

The radio-opaque dots to the right of the inlet connector are used to identify the pressure range of the valve:

		
SM1-L (Low Pressure)	SM1-M (Medium Pressure)	SM1-H (High Pressure)
1 radio-opaque dot ●	2 radio-opaque dots ●●	3 radio-opaque dots ●●●

7.4. PATENCY TEST (POST-OPERATIVE)

There are two steps for the post-operative test on the patency of the shunt :

7.4.1. Testing the patency of the ventricular catheter

NOTE

This control is possible with SM1A (antechamber) and SM1B (burr hole reservoir) type models. For the SM1 type models, a ventricular catheter with reservoir must be used.

Pinch the catheter with a finger just after the valve outlet connector.

With another finger, press the reservoir to make the CSF flow back into the ventricular catheter. A reservoir that cannot be compressed easily or does not fill quickly may indicate there is an obstruction in the ventricular catheter.

7.4.2. Patency test downstream of the reservoir (valve and distal catheter)

NOTE

This control is impossible with a SM1B type model (burr hole reservoir) because there is no access to the ventricular catheter upstream of the reservoir.

Pinch the catheter with a finger just before the reservoir, then with another finger press the reservoir to push the CSF through the valve and distal catheter. A reservoir that cannot be compressed easily may indicate an obstruction either of the valve or the distal catheter.

NOTE

Do not rely only on the characteristics of the patency test to diagnose an obstruction in the shunt system. Obstruction of a shunt system can occur in any of its components and should be diagnosed first of all by the clinical data and additional examinations.

7.5. SAMPLING THE CSF AND INJECTION

Access to the CSF is obtained by pricking the reservoir with a 24G (or smaller diameter) Huber needle.

The integral reservoir on the SM1A and SM1B models is designed for occasional use.

Its watertight performance is reduced after very frequent pricking into the dome.

- To inject in the proximal direction, compress the catheter just after the valve outlet connector.
- To inject in the distal direction, compress the catheter upstream of the reservoir.

NOTE

Elective injection in the distal direction is not possible with a model of the SM1B type (burr hole reservoir) because there is no access to the ventricular catheter upstream of the reservoir.

PRECAUTIONS

Prior to the injection enquire about the compatibility of the injected product with silicone.

Ensure that the base of the reservoir is not crossed with the needle.

Do not inject into, or take samples from, the CSF without having tested the shunt patency. Significant overpressure could damage the shunt if it is known that there is an obstruction.

Ensure that the orifice of the needle is inside the reservoir.

Do not inject too fast or inject too great a volume. The increase in pressure could damage the shunt.

Do not use a syringe with a volume of less than 10cc for injections or taking samples. Too great a pressure could damage the shunt.

8. Precautions for the Daily Life of the Patient

A Patient Identification Card (PC-SM1) is supplied with the Sophy® Mini valve. It enables the neurosurgeon to consult and update information relating to the implanted device (reference, operating pressure, implantation site, etc.) systematically and to ensure that the illness is properly monitored.

WARNING

The doctor is responsible for informing the patient or his/her family that the person fitted with a shunt must avoid any activity that may subject this shunt to direct shocks (violent sports, etc.) as these are likely to damage it.

PRECAUTION

The patient must be warned that vibrations due to the CSF flow may possibly be felt because of the implantation of the valve on the skull.

9. Complications / Side effects

Complications which may result from the implantation of a CSF shunt system include the inherent risks in the use

of drugs, any surgical intervention and the insertion of a foreign body.

NOTE

Patients treated with a shunt system must be closely monitored post-operatively in order to detect any signs of complications early.

The doctor is responsible for educating the patient or his/her family about CSF shunt systems, in particular describing the complications linked to implanted shunt systems as well as giving explanations about possible alternative therapies.

The main complications of shunts are obstruction, infection and over-drainage. These complications require the rapid intervention of a doctor.

9.1. OBSTRUCTION

Obstruction is the most frequent complication in shunt systems. It can occur at any point in a shunt.

The ventricular catheter can be obstructed by a blood clot, cerebral tissue or even tumoral cells.

The end of the ventricular catheter can also become embedded in the choroid plexus or in the ventricular wall, either directly or following a collapse of the walls, a consequence of over-drainage.

The cardiac catheter can be colonized by a thrombus while the appearance of a clot around the catheter could cause an embolism in the pulmonary circulation.

The peritoneal catheter may become obstructed by the peritoneum or by intestinal loops.

Loss of patency in a shunt may also be the result of an obstruction by fragments of cerebral tissue or by biological deposits (protein deposits, etc.).

Obstruction of the shunt will quickly result in the reappearance of the signs and symptoms of intracranial hypertension.

These signs and symptoms vary from patient to patient and over time.

In infants and young children, the symptoms may be an abnormal increase in the size of the skull, a bulge in the fontanelles, dilation of the scalp veins, vomiting, irritability with a lack of attention, downward deviation of the eyes, and sometimes convulsions.

In older children and adults, intracranial hypertension due to hydrocephalus may be the cause of headaches, vomiting, blurred vision, diplopia, drowsiness, slowing of movements, gait disorders or psychomotor slowing which could lead to total invalidity.

If an obstruction is confirmed and a patency test does not make it possible to reduce the obstruction, revision surgery or removal of the device must be envisaged.

9.2. INFECTION

Chronic malfunction of the shunt could cause a leak and a discharge of CSF along its length increasing the risk of infection.

Local or systemic infection is another possible complication of CSF shunt systems. It is generally

secondary to the colonization of the shunt by cutaneous germs. Nevertheless, as for all foreign bodies, any local or systemic infection can colonize the shunt. Erythema, edema and skin erosions along the length of the shunt may be an indication of an infection of the shunt system.

Prolonged, unexplained fever may also be the result of a shunt system infection.

Septicemia, favored by an alteration in general status, can start from a shunt infection.

If there is infection, removal of the system is indicated in conjunction with the start of a specific treatment by a general or intrathecal route.

9.3. OVERDRAINAGE

Overdrainage can result in a collapse of the ventricles (slit ventricle syndrome) and the appearance of a subdural hematoma.

In children, depression of the fontanelles, overlapping of the scalp bones, even a craniostenosis or a change from communicating hydrocephalus to obstructive hydrocephalus by stenosis of the Aqueduct of Sylvius could occur.

Adults can present with a variety of symptoms such as vomiting, auditory or visual disorders, drowsiness or even headaches in the upright position but which improve in the supine position.

However, immediate drainage of a subdural hematoma may be indicated.

9.4. OTHER

Failure of a shunt system may also be linked to disconnection of its various components.

The ventricular catheter may migrate inside the ventricle. The peritoneal catheter may migrate into the peritoneal cavity under the action of the peristaltic waves of the intestine, while an atrial catheter may migrate into the right-hand cavities of the heart following the blood flow.

Perforation or occlusion of abdominal viscera by the peritoneal catheter could occur.

Body growth may progressively cause the catheters to exit their atrial or peritoneal insertion sites.

These malfunctions require the shunt to be repositioned immediately.

Cases of cutaneous necrosis over the implantation site are possible.

Over time fibrous adhesions may fix the ventricular catheter in the choroid plexus or the cerebral tissue. If removal is being considered, gentle rotation of the catheter about its axis may make it possible to free it. The catheter should never be withdrawn forcibly. If it cannot be taken out without forcing, it is preferable to leave it in place rather than risk an intra-ventricular hemorrhage.

Cases of allergy to silicone have been described.

Cases of epilepsy after implantation of a ventricular shunt have been described.

The ruby ball in the valve can potentially take up an off-center position on its housing due to the presence of a

cluster of cells or protein deposits. Among others, such situations can cause :

- loss of regulatory function in the valve potentially increasing the risk of overdrainage.
- an impaired anti-reflux function.

Blood clots, cerebral cells or tumoral cells contained in the CSF could lodge in the valve mechanism, which would have the potential to cause changes in the operating characteristics of the valve.

10. Storage conditions

Keep in the original packaging.

Keep in a cool, dry place, away from sun light and heat.

11. Processing of the products after use

11.1. DESTRUCTION AFTER USE

An unpacked, used or explanted Sophy® Mini monopressure valve must be destroyed in accordance with the procedures in force in the health establishment.

11.2. RETURN OF PRODUCTS

If an explanted valve needs to be returned to Sophysa for analysis, it must be returned immersed in water, indicating if necessary whether cleaning has been performed.

Never use a saline solution likely to form deposits in the valve body which could block the rotor.

In order to assess the returned product properly, it must be accompanied by an explanatory Return to Manufacturer Authorisation form.

12. Warranty

Sophysa guarantees that this medical device is free of any material and manufacturing defects. Apart from this guarantee, Sophysa does not provide any other guarantee, express or implicit, including commercialization or adaptation for a particular use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice occurring directly or indirectly from the use of this device. Sophysa does not authorize anyone whomsoever to take responsibility on its behalf for its products.

The performance of Sophy® Mini monopressure valves is only guaranteed with the range of silicone catheters and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of catheters provided that their internal diameter is identical to that of the catheters recommended by Sophysa.

13. Symbols

REF

Catalog reference

Manufacturer

Date of manufacture

Caution, see the Instructions for Use

STERILE EO

Sterilized using Ethylene Oxide

Do not re-use

Do not re-sterilize

Use by

LOT

Batch code

SN

Serial number

CE Conformity Marking

14. References

14.1. SOPHY® MINI MONOPRESSURE VALVES FOR CSF SHUNTING

SM1-L	Sophy® Mini monopressure valve, Low Pressure Fixed pressure valve, 50 mmH ₂ O
SM1-M	Sophy® Mini monopressure valve, Medium Pressure Fixed pressure valve, 110 mmH ₂ O
SM1-H	Sophy® Mini monopressure valve, High Pressure Fixed pressure valve, 170 mmH ₂ O
SM1A-L	Sophy® Mini monopressure valve, Low Pressure, with antechamber Fixed pressure valve, 50 mmH ₂ O. Integral antechamber.
SM1A-M	Sophy® Mini monopressure valve, Medium Pressure, with antechamber Fixed pressure valve, 110 mmH ₂ O. Integral antechamber.
SM1A-H	Sophy® Mini monopressure valve, High Pressure, with antechamber Fixed pressure valve, 170 mmH ₂ O. Integral antechamber.

SM1B-L	<p>Sophy® Mini monopressure valve, Low Pressure, with “burr hole” type reservoir</p> <p>Fixed pressure valve, 50 mmH₂O. Integral “burr hole” type reservoir.</p>
SM1B-M	<p>Sophy® Mini monopressure valve, Medium Pressure, with “burr hole” type reservoir</p> <p>Fixed pressure valve, 110 mmH₂O. Integral “burr hole” type reservoir.</p>
SM1B-H	<p>Sophy® Mini monopressure valve, High Pressure, with “burr hole” type reservoir</p> <p>Fixed pressure valve, 170 mmH₂O. Integral “burr hole” type reservoir.</p>

Availability may vary according to country.

14.2. SHUNT KITS FOR CRANIAL IMPLANTATION

SM1-2010L	<p>Complete Sophy® Mini monopressure kit, Low Pressure</p> <p>Sophy® Mini monopressure valve, Low Pressure (50 mmH₂O) with preconnected peritoneal catheter (B905S). Right ventricular catheter (BO19-10).</p>
SM1-2010M	<p>Complete Sophy® Mini monopressure kit, Medium Pressure</p> <p>Sophy® Mini monopressure valve, Medium Pressure (110 mmH₂O) with preconnected peritoneal catheter (B905S). Right ventricular catheter (BO19-10).</p>
SM1-2010H	<p>Complete Sophy® Mini monopressure kit, High Pressure</p> <p>Sophy® Mini monopressure valve, High Pressure (170 mmH₂O) with preconnected peritoneal catheter (B905S). Right ventricular catheter (BO19-10).</p>
SM1A-2010L	<p>Complete Sophy® Mini monopressure kit, Low Pressure, with antechamber</p> <p>Sophy® Mini monopressure valve, Low Pressure (50 mmH₂O), integral antechamber, with preconnected peritoneal catheter (B905S). Right ventricular catheter (BO19-10).</p>
SM1A-2010M	<p>Complete Sophy® Mini monopressure kit, Medium Pressure, with antechamber</p> <p>Sophy® Mini monopressure valve, Medium Pressure (110 mmH₂O), integral antechamber, with preconnected peritoneal catheter (B905S). Right ventricular catheter (BO19-10).</p>
SM1A-2010H	<p>Complete Sophy® Mini monopressure kit, High Pressure, with antechamber</p> <p>Sophy® Mini monopressure valve, High Pressure (170 mmH₂O), integral antechamber, with preconnected peritoneal catheter (B905S). Right ventricular catheter (BO19-10).</p>

Technical specifications and List of product references may be modified without notice. For all product dimensions indicated in this document, tolerances are: ± 5 %.

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