


PATIENT POSITIONING FOR THE RADIOGRAPHIC EXAM AND LANDMARKS OF THE VALVE



Sophy® is a trademark registered by Sophysa


©Sophysa 2021

DFSh-SM8-02(EN) Rev000 - 01/2021



SOPHY® Mini SM8

READING THE OPERATING PRESSURE



P1

P2

P3

P4

P5

P6

P7

P8

SM8-140

SM8

SM8-300

SM8-400

10 mm H₂O

30 mm H₂O

50 mm H₂O

80 mm H₂O

25 mm H₂O

50 mm H₂O

75 mm H₂O

120 mm H₂O

40 mm H₂O

70 mm H₂O

100 mm H₂O

150 mm H₂O

60 mm H₂O

90 mm H₂O

125 mm H₂O

190 mm H₂O

80 mm H₂O

110 mm H₂O

150 mm H₂O

230 mm H₂O

100 mm H₂O

140 mm H₂O

180 mm H₂O

270 mm H₂O

120 mm H₂O

170 mm H₂O

220 mm H₂O

330 mm H₂O

140 mm H₂O

200 mm H₂O

300 mm H₂O

400 mm H₂O